

Panasonic

REAL PRO **ULTRA**TM
Prestige

EPMAJ7K



OUR DNA

Panasonic

Introducing the **Real Pro Ultra Prestige Massage Experience by Panasonic**. The perfect marriage of luxury design and powerful, yet intuitive massage techniques. The **EPMAJ7K** is designed to help you achieve therapeutic muscle relaxation and increased blood flow throughout your body to promote good health. Our newly developed Lower Back & Glute massage program was created to specifically target and help alleviate the effects of lower back pain.

Our unique 3D Massage System is comprised of 3 collaborative elements; Thermal Massage Technology, Professional Massage Techniques and Professional Stretching Techniques, which work in concert to produce a professional quality massage experience. Experience the freedom of daily personalized massage therapy in the comfort of your own home.



Experts in massage therapy since 1936
80 Years of Development & Discovery

MASSAGE ELEMENTS

3D Massage Elements


Thermal
Massage
Technology



Professional
Stretching
Techniques



Professional
Massage
Techniques



Enjoy a professional massage
in the comfort of your own home.

Panasonic's latest 3D Massage System simulates the
warm, soothing hands of a professional masseuse.

Our Unique 3D Massage System

uses a high performance brush-less motor and allows the rollers to move up and down, left to right and front to back.

Our Automatic Intensity Sensor

intelligently adjusts applied pressure for maximum comfort.

Our Intuitive Micro Processor

allows for excellent responsiveness to replicate intricate Shiatsu techniques.

Our High Speed Digital Control

accurately replicates multiple professional massage techniques.



JUNETSU HEATED MASSAGE

The healing power of heat.
Relax, and revive your body with a heated massage.

The calming strokes of a masseur’s warm, skillful hands are replicated by our high performance, built-in thermal massage rollers. These skillful movements help ease swelling and relieve overall tension.

Junetsu Massage
Panasonic’s Exclusive Junetsu Technique specifically targets acupoints found in the back and neck with a deeply invigorating circular motion designed to energize tired muscles.

Shiatsu Heated Body Massage
Our unique 3D massage system replicates intricate Shiatsu Massage techniques by applying heat and pressure to targeted acupuncture points to promote total body relaxation.

NEW! Foot Roller Massage Technology
Massage rollers activate the arch of the foot to alleviate stiffness and use the soothing pressure of nodules to loosen the ball of the foot, while alternating airbags apply gentle pressure for maximum relaxation.



NEW! LOWER BACK & GLUTES

NEW! Lower Back & Glute Massage

Targeted therapy to help relieve lower back pain, stiffness and tension

The lower lumbar and glute muscles work together to support the spine and lower back. The **EPMAJ7K** features new targeted massage programs specifically designed to relieve tension in the lumbar and glute regions, to help provide relief from persistent back pain.

Choose from a wide range of customized options ready to address your specific needs. For a more intense massage experience, raise the ottoman up to a comfortable position.



NEW! SHIATSU FOOT ROLLERS

From your head to your feet
Relieve Tired Feet with our **NEW** Luxury Foot Massage Program



NEW!

Foot Roller Massage Technology


The **EPMAJ7K's** innovative foot massage technology activates the arch of the foot with a firm rolling motion, smoothing away tension and fatigue and bringing relief to sore arches.

Shiatsu Foot Massage

The targeted pressure of massage nodules help activate acupoints in the feet to loosen the ball of the foot, while alternating airbags apply gentle pressure for maximum relaxation.

Invigorating Airbag Massage

The **EPMAJ7K** massaging ottoman features both calf and foot airbag systems to help enhance the massage experience by temporarily increasing blood flow for an energizing feeling.



NEW! SPIRAL KNEAD MASSAGE


NEW! Spiral Knead Massage

**Innovative Spiral Kneading Massage Technique,
a Panasonic EPMAJ7K Exclusive Feature**

The **EPMAJ7K's** state-of-the-art thumb-like massage rollers relieve and loosen tight muscles and knots by executing a new spiral kneading technique delivering an intense yet comfortable graduated circular motion that invigorates tired muscles.

The Personal Touch

You can design your own personal luxury massage experience with 6 different types of spiral knead programs:



STRETCH MODE

NEW! Professional Stretching Techniques

Improve Flexibility, Circulation and Reduce Tension

Our combined 3D massage techniques and full-body airbags stretch deep into the muscle tissue for a powerful, therapeutic stretching massage **with a new focus on the lower back and core.**

Neck, Chest, Pelvis & Leg Stretch Techniques

Airbag technology elevates and secures these core areas while rollers gently massage. As the airbags slowly deflate, they create a relaxing stretch along with intermittent kneading for a pleasurable experience.

NEW! Lower Back Stretch

Airbags inflate to secure and lift the lower body while massage rollers gently push up and stretch the lower back. Air pressure is then reduced in the seated airbags, stretching the surrounding area.

NEW! Core Stretch

This new feature stretches the muscles of the abdomen, lower back and glutes. These major muscle groups work together to support healthy posture, with the application of heat and stretching – the body relaxes from the core.



ADDITIONAL FEATURES

Full Body Air Bag Massage
Temporarily Stimulates Blood Flow for Invigorating Energy

Features

Armchair Use	Relax in a comfortable, upright, seated position
Flat Recline 170°	170° of control providing custom personal comfort
Premium Synthetic Leather	Attractive cross-hatch design is easy-to-clean
Wheels	Moving made easy
Safety Lock	Power and lock switch prevent malfunctions
Height Adjustable	Ottoman extends up to 18 cm



AUTOMATIC MASSAGE PROGRAMS

Design Your Own Personalized Massage Experience



More intelligent control with 3 memory settings for the whole family. Our **NEW! Massage Strength Control** allows you to intensify or soften the standard massage strength setting for all automatic programs.

Full Body Programs		
PROGRAM	INTENSITY	DURATION
SHIATSU	STRONG	19 MIN
DEEP	STRONG	19 MIN
STRETCH	STANDARD	16 MIN
SWEDISH	SOFT	19 MIN

Targeted Programs		
PROGRAM	INTENSITY	DURATION
NECK / SHOULDER	STANDARD	16 MIN
UPPER BACK	STANDARD	16 MIN
CORE	STANDARD	16 MIN
LOWER BACK / HIP	STANDARD	16 MIN
QUICK	STRONG	10 MIN


MANUAL MASSAGE PROGRAMS

Manual Massage Programs

Customize your massage with over 100 different combinations


Area	Ultra Knead	Shiatsu	Knead	Swedish	Rolling	Spiral
Neck						
Shoulder						
Middle Back						
Lower Back						
Full Back						
Glutes						

Wider Neck & Shoulder Range



Massage rollers extend up to 10 cm, ensuring effective relief of shoulder stiffness. A wide area of therapeutic massage is covered, up to 79 cm vertically and up to 25 cm horizontally.

Auto Body Scanning System



The sophisticated auto body scanning system is able to detect the prominent points of your body and adjust the massage rollers to fit your unique structure. The sensor technology delivers a customized massage specifically tailored for each user.

EPMAJ7K SPECIFICATIONS

Power Supply		AC 120 V 60 Hz
Rated Current		3.0 A
Upper Body Massage	Massage Area (up-down)	Approx. 79cm (The massage heads move a total of approx. 63 cm)
	Massage Area (left-right)	Distance between massage heads during operation (including width of massage heads) Neck, shoulder, back: Approx. 7 - 21 cm Lower back & hip: Approx. 7 - 25 cm Distance between massage heads (including width of massage heads) during back roll and up-and-down movement: Approx. 7 - 18 cm
	Massage Area (forward-back)	Projected range of massage head travel (intensity adjustment): Approx. 10 cm
	Massage Speed	Knead: From approx. 4 - 73 times/min Tap: From approx. 190 - 520 times/min Moving Speed: From approx. 1 - 5 cm/sec
Air Massage	Air Pressure (Calf)	Level 1: Approx 17 kPa, Level 2: Approx. 32 kPa, Level 3: Approx. 35 kPa (May vary with the area slightly.)
Sole Massage	Massage Speed	From approx. 8 - 21 times/min
Reclining Angle		Backrest: Approx. 120° - 170° Leg rest (leg & sole massage section): Approx. 5° - 85°
Timer		Up to 19 minutes (automatically shuts off)
Dimensions	When not reclined	(H x W x D) Approx. 115 x 90 x 130 cm
	When reclined	(H x W x D) Approx. 72 x 90 x 200 cm
Weight		Approx. 87 kg
Material		High Quality Synthetic Leather
Preset Programs		Swedish, Deep, Shiatsu, Stretch, Quick, Core, Neck/Shoulder, Upper Back, Lower Back/Hip